


# • THE SUPERPOLE SYSTEM



**A**

## **A SUPERPOLE BASIC**

**Ease of Installation** - Easy to install or to relocate, no drilling or screws, no marring of floors/ceilings.

**Flexible** - Install anywhere where there are floors or ceilings. Place support where it is needed.

**Flexible** - Custom heights of 120"/140" ceilings (see also Uni-fit extender).

**Attractive** - Clean and modern design

**Solid** - Made of heavy gauge steel for rigid support.

**Comfort** - High-quality grip, especially for those with limited grip strength

**Secure** - Rubber pad interface top and bottom

**Toe Friendly** - Small and low profile base plate.


## **B SUPERBAR Pivoting and Locking Rail**

**Safety** - Pivots freely to move with you, step by step (always nearby for safety).

**Safety** - Can lock every 45° to assist you with standing or sitting.

**Ease of Use** - Unlocks easily (lift 0.25" to move freely).

**Safety** - Bar locks automatically every 45° (i.e. in the event of loss of balance).



**B**

**Safety** - Can prevent a forward fall by locking (common cause of broken hip).

**Safety** - Offers a wider coverage range due to 16" wide pivoting handle.

**Safety** - Full support for standing and 2-3 step transfer (every step of the way).

**Safety** - More stable support due to wider support stance.

**Safety** - Easier weight bearing on the horizontal surface of the horizontal hand rail.


**Ease of Use** - Less range of motion required by client (no more reaching for support).

**Flexible** - Fits in small areas such as bathrooms - more compact than walker.

**Versatile** - SuperBar can pivot into narrow passages such as hallways.

**Bed Application** - Allows client to weight shift, sit, pivot, stand and transfer.

**Bathtub Application** - SuperBar supports step-over bathtub wall while reaching for wall grab bars. With glass doors, where clamp-on type hand rail cannot be installed, SuperBar can pivot in/out to allow doors to close while showering.



**C**

**Toilet Application** - For pulling up safely to standing position. Often will service toilet and adjacent bathtub.

## **G SUPERTRAPEZE with Easy-Up! Handle**

**Ease of Use** - combination of trapeze & bedside pole. Only product to offer complete support for sitting up/pivoting/standing up from bed.

**Effective** - Includes EasyUp! trapeze handle with ladder-like set of 2 rails, allowing client to more easily "climb up" to sitting position.

**Flexible** - Support arm adjustable in height / angle. Strap adjustable in length to achieve best height and location of the trapeze.


**Ease of Use** - Open design of trapeze handle allows it to be hooked onto pole for storage

**Ease of Use** - Does not need special bed or mounting brackets - works with any bed.

**Versatile** - Does not force bed outward from wall as with conventional trapeze systems.

**Value** - More functions for less money.

**Versatile** - Will install beside water and electrical beds unlike floor mounted trapeze.



**D**

## **D SUPERTRAY**

**Ease of Use** - Quickly installs on or removes from SuperPole.

**Versatile** - Provides function of an over-bed table, moves easily and does not occupy floor space.

**Versatile** - Great for meals, reading material, medications etc.

**Ease of Use** - Slim profile allows for easy cleaning and storage.

**Ease of Use** - Pivots and removes from pole. Tray removes from tray frame.

Visit our showroom at the BARRIER FREE LIVING CENTRE, 2700 Dufferin St. #24, Toronto ON Canada M6B 4J3



**BARRIER FREE**  
ARCHITECTURALS INC.



Architectural Products for Barrier Free Living

Tel 416 783 5331 Fax 416 783 6505

Toll free 1 877 717 7027 e-mail info@barrierfree.org

web site  
[www.barrierfree.org](http://www.barrierfree.org)